



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LEISEGANG 1 Colposcope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Leisegang 1 is a chair-mounted colposcope for gynaecological examinations,
- Mounted on a 2-part arm,
- Light source: 6 V / 20 W halogen bulb,
- Features a mechanism that allows the colposcope to be comfortably adjusted in height within a range of 33 cm,
- Equipped with a handle for better control of the device,
- Infinitely variable manual adjustment:
 - The height of the colposcope head on the pole within 15 cm,
 - Light intensity,
 - Pupil distance and eyepiece correction,
- Dimensions: 90 x 29 x 146 cm,
- Weight: 20 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

