



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER VAPOR 19.3 Vaporiser for anaesthetic agents (Halothane) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Germany,
- The Vapor 19.3 is an anaesthetic vaporiser, it allows efficient and economical administration of the anesthetic agent,
- Product features:
 - Designed for Halothane,
 - Vaporiser volume: max 200 ml,
 - Operating temperature: 15° C - 35° C,
 - Operating pressure: 700 to 1013 mbar,
 - Infinitely adjustable anesthetic concentration,
 - Unique transport function - vaporisers can be carried in any position. The transport safety device prevents spillage of the anaesthetic so that it can be used immediately after transport,
 - All filling systems are specifically designed for specific inhalation anaesthetics. This means that the risk of mistakenly using the wrong anaesthetic is eliminated,
 - The Dräger attachment (DW-2000) is used in this vaporiser,
 - Fill type: Funell Fill,
- Included:
 - User manual in English (PDF),
 - Overall dimensions: 16.5 x 10 x 25 cm,
- Weight: 7 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

