



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ESCHMANN VP-35 Medical suction unit (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in England,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 0.7 A,
- VP35 is a high vacuum high flow suction device,
- Specifications:
 - High flow rate of 35 liters/min,
 - High maximum vacuum of 600 mm Hg (- 80 kPa),
- Equipped with:
 - Built-in pressure gauge to adjust suction power,
 - Antibacterial, hydrophobic filter,
 - Convenient handle for transporting the suction unit,
 - Mobile on four casters,
- Includes:
 - Two brand new 2000 ml autoclavable jars from SEERES,
 - Two brand new 2000 ml disposable cartridges,
 - Suction hose,
- Overall dimensions: 41 x 44 x 67 cm,
- Weight: 14.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

