



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **VERMEIREN V300 30° COMFORT Wheelchair (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Belgium, (2018),
- V300 30° Comfort is a wheelchair with head and back stabilization made of high-quality aluminum, with advanced adjustment options,
- Advantages:
  - Specialized Contour type upholstery offers comfort to the sitting person at the highest level,
  - Deep side contouring of the backrest perfectly stabilizes the torso,
  - Seat contours protect the pelvis from tilting,
  - ORTHO type headrest perfectly supports and stabilizes the head,
  - Adjustable footrest with calf support ORTHO type - is characterized by the fact that the sitting person can not only set a different angle, but also can accurately adjust the positioning of individual supports, which are directly responsible for the correct posture of the person sitting in the wheelchair as well as her/his correct stabilization:
    - Calf support adjustable in depth and height,
    - Foot support adjustable in height as well as adjustable foot support angle,
  - This wheelchair is dedicated for people with unstable torso, especially for patients with paralysis or quadriplegia, triple limb paresis, post-stroke and MS patients,
- Product characteristics:
  - Rear wheels detachable with quick-release couplers,
  - Rear wheels in an adapter to change the center of gravity,
  - Contoured PU SOFT upholstery - L34L35,
  - Special hooks attaching the backrest to the frame of the wheelchair allow the sitting person to move the backrest vertically and horizontally. In this way the user can raise the backrest or increase the depth of the seat,
  - 3D orthopedic headrest - L58,
  - Adjustable footrests with calf support -BZ7 ORTHO,
  - Height-adjustable armrests - B05,
  - Variable angle footrest plate,
  - 3 seat heights from the ground,
  - 3 heights of elbow support,
  - 3 depths of elbow support,
  - 5 levels of seat tilt,
  - Variable angle fork,
- Folding and transporting the wheelchair is very easy, it can be folded in 2 two minutes after removing the upholstery and wheels. Thanks to a unique design solution, the wheelchair folds to a very minimal size and can be transported even in a small trunk,
- Included:
  - User manual in Polish (PDF),
- Dimensions:
  - Seat: 50 x 42 cm,
  - Backrest (without headrest): 43 x 53 cm,
  - Overall: 104 x 62 x 117 cm,
- Maximum user weight: 130 kg,
- Weight: 30 kg,
- Cat 01,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

