



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS OEV191H LCD Medical monitor (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 125 VA,
- Olympus OEV191H is an HDTV quality compliant medical monitor,
- Progressive scanning reduces image distortion and flicker,
- Anti-glare coating minimizes light reflection,
- Image enhancement function,
- 19" TFT LCD matrix,
- Horizontal and vertical viewing angle: 85°,
- Resolution: 1280 x 1024 px,
- Menu in English,
- Keypad on control panel in language: English,
- Features the following inputs and outputs:
 - COMPUTER,
 - VIDEO IN-OUT, Y/C IN-OUT,
 - REMOTE: REMOTE 1, RS-232C,
 - RGB/COMPONENT: A (G/Y IN-OUT, B/PB IN-OUT, R/PR IN-OUT, EXT SYNC IN-OUT), B (G/Y IN, B/PB IN, R/PR IN, EXT SYNC IN),
- Dimensions: 43 x 38.5 x 10 cm,
- Weight: 8 kg,
- Cat 03,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински



