



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EURONDA Model E9 MED Class B autoclave (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Italy,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1400 W,
- To guarantee the highest quality of the sterilization process, the Euronada E9 MED autoclave is equipped with a steam generator. This produces highly saturated steam in a short period of time. The steam passage to the boiler is protected by a filter, which can be easily removed for cleaning. The cover of the water tanks can be easily removed and the clean and used water tanks can be cleaned. To extend the life of the vacuum pump, a special separator has been installed to eliminate water getting into the pump,
- The memory allows recording the last 40 cycles, the results of which can be read (printed) in order or by selecting, for example, those that were successful or those that were aborted,
- The autoclave can be equipped with an Aquafilter water conditioner (not included), dispensing demineralized water automatically. The conditioner has a water quality sensor, extending the life of the autoclave. To improve handling, the used water tank can be connected to the sewage system,
- E9 MED is equipped with an electro-magnetic door lock to prevent accidental opening of the chamber door.
- The boiler is made of high-quality stainless steel,
- The rack inside the chamber can accommodate four trays or, when inverted, three cassettes,
- The control panel of the autoclave is equipped with a large liquid-crystal display. With only 3 buttons, the operator can easily select the desired functions and settings of the autoclave. An additional advantage is the menu in Polish language,
- Features a "Process Controller" function, which informs the user in real time about the process: temperature, pressure and time during all phases of the selected program. The parameter reading is converted into a message to the user with regard to: safety functions, batch conditions inside the chamber, maintenance information, autoclave memory reading, printing functions and functions of external devices such as water conditioner,
- Test cycles: Vacum Test, Helix Test, B&D Test,
- The autoclave itself intelligently prompts on the display when to: clean the boiler, replace the door seal, bacteriological filter, when to clean the water tanks or when to perform a technical inspection - saving on service,
- Very quiet, only 50 dB - important when the autoclave is in the office,
- Capacity: 18 liters,
- Built-in thermal printer,
- Additional archiving devices can be connected to the autoclave (3,000 E-MEMORY cycles) - not included,
- Selectable menu in language: PL, IT, DE, EN, ES, FR, LT, PT, TR,
- Included:
 - User manual in English (PDF),
 - Basket, two trays and handle,
 - Drainage tube,
 - Power cord,
- Dimensions: 61 x 45 x 44 cm,
- Weight: 41 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

