



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DEL TOM WPC-220 BMH nursing and rehabilitation bathtub (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- DELTOM WPC-220 BMH bathtub for hydrotherapy and rehabilitation. It is a free-standing device that does not require housing and foundation which allows it to be set up freely without the need for special room adaptation,
- Advantages:
 - Equipped with hydromassage that gives the opportunity to relax,
 - Accessible from all sides and adjustable height of the bathtub allows medical personnel to assume a comfortable, ergonomic position,
- Equipped with:
 - Mixable water faucet,
 - Shower,
 - Digital thermometer for measuring water temperature,
 - Underwater hydromassage system,
 - Adjustable hydromassage jet - its direction and intensity,
 - Underwater spotlight illuminating the bathtub,
 - Seat with headrest,
 - Connection to cold and hot water,
- Height adjustment of the bathtub by wired remote control at the level of the edge of the basin: 79 - 114 cm,
- Maximum capacity of the basin: 200 L,
- Dimensions:
 - Overall: 100 x 200 x 120 cm,
 - Basin: 60 x 120 x 47 cm,
- Weight: 180 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

