



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG BAEHR Vento 75 Dry heat sterilizer (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 650 W,
- MELAG BAEHR Vento 75 is a hot air sterilizer, fully safe and durable device. It meets all the requirements for compact devices,
- Safety is guaranteed by the function of automatic shutdown of the sterilization cycle including turning off the power after the cycle is completed,
- The device can be used in veterinary and medical fields,
- Dry sterilization,
- Forced air circulation,
- Capacity: 4.3 L,
- Infinitely adjustable temperature from 60° C to max. 180° C,
- Features a timer switch adjustable from: 0 - 120 min,
- Stainless steel housing, strong construction and long life of the device,
- The set includes:
 - Tray,
 - Tray holder,
 - User manual in German (PDF),
 - Power cord,
- Dimensions:
 - Tray: 18 x 28 x 2 cm,
 - Chamber: 19.5 x 29 x 7.3 cm,
 - Overall: 31 x 41.5 x 26 cm,
- Maximum load of tray: 2 kg,
- Weight: 12 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

