



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## ZALMED AAD-80 Diagnostic audiometer (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 100 VA,
- Automatic and manual hearing test function,
- Two channels,
- Modulation frequency: from 0.2 to 4 Hz,
- Modulation depth from 0.25 to 4 Hz,
- Frequency from 250 Hz to 8000 Hz,
- SISI function,
- Built-in speaker,
- Control panel keyboard in Polish,
- Included:
  - Power cord,
  - Tonsil SD-307 headphones,
  - B-71 bone conduction headphone,
  - Patient buttons,
- Dimensions: 47 x 47 x 21 cm,
- Weight: 16 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински



