



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STRYKER Arthroscopy pump (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power consumption: 70W,
- STRYKER arthroscopy pump, with its simple and intuitive operation, improves the comfort and safety of work, and thus provides full functionality,
- Maximum pumping pressure: 500 mmHg,
- Flow range: 0 - 1000 ml / min,
- Included:
 - Power cord,
- Dimensions: 36 x 15 x 34 cm,
- Weight: 9.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

