



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MED X CHANGE EVOLUTION4K EVO4K200 4K video recorder (Reconditioned)

- New (demo),
- Technical condition: perfect,
- Visual condition: perfect,
- Made in USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Med X Change's Evolution4K EVO4K200 offers 4K UHD video recording in a compact unit, while providing superior 4K UHD video quality designed for both Apple and Windows editing,
- The Evolution4K allows users to enter patient data, capture 4K UHD video and 4K images, play back 4K video, automatically print to medical printers, and record to external hard drives or USB 3.0 flash drives. The device is compatible with medical foot switches (RS-232), 4K surgical microscope and 4K endoscopic cameras (2 x 3.5mm triggers), encrypts all ePHI data at rest and in transit, and can store data on a large internal hard drive,
- Operated via an integrated 5-inch capacitive multi-touch screen that uses swipe, zoom gestures, similar to the latest mobile devices,
- Specifications:
 - Full 4K UHD I HD video/photo recording:
 - 3G-SDI model (4, EVO4K200),
 - 4K UHD 60 Hz input: 3G-SDI,
 - 4K UHD 60 Hz output: 3G-SDI,
 - 4096 x 2160p or 3840 x 2160p 50/60 Hz,
 - Recorded/captured 4K UHD files: H.264, .mp4, JPEG, PNG,
 - Integrated 12.7 cm LCD touch screen,
 - Control all software from integrated touch screen,
 - Multi-touch / swipe interface,
 - Duplicate the internal diagonal touchscreen on any HID touchscreen monitor via VGA connector,
 - Image capture:
 - 4K: 4096 pixels x 2160 pixels = 9 MP,
 - UHD: 3840 pixels x 2160 pixels = 8 MP,
 - HD: 1920 pixels x 1080 pixels = 2 MP,
 - Typical supported camera/video resolutions:
 - 4K and 4K UHD: 4096 x 2160p and 3840 x 2160p,
 - HD: 1080p, 1080i, 720p,
 - 4K and 4K UHD Hz: 60, 59,94, 50,
 - 1080p Hz: 60, 59,94, 50,
 - 1080i Hz: 60, 59,94, 50,
 - 720p Hz: 60, 59,94, 50,
 - 3D HD Side by Side: 1080p, 1080i, 720p,
 - Top/bottom 3D HD: 1080p, 720p,
 - Patient input,
 - One-touch recording,
 - Simultaneous recording to internal hard drive and 2 USB 3.0 locations,
 - DICOM integration,
 - Live video transmission over IP via RTSP protocol,
 - Printing, playback, copying, viewing of recorded material,
 - Built-in 1 TB hard drive,
 - Linux operating system,
- Selectable menu in EN, DE, ES, FR, IT, NL, PR, FI, TR, CZ, RU, JP, KR,
- Included:
 - Power cable,
 - 6x SDI 4K - BNC connectors,
- Dimensions: 28 x 21 x 11 cm,
- Weight: 6.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN DUTCH
Мы говорим по русски
Мы говорим по украински

