



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EPENDORF CENTRIFUGE 5702 R Laboratory centrifuge (Ex-Display)

- New, exc-display,
- Technical condition: perfect,
- Visual condition: perfect,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power consumption: 380 W,
- The Centrifuge 5702 R is a refrigerated centrifuge model for heat-sensitive samples,
- Built-in backlit LCD display,
- Keypad on the control panel in English,
- Product features:
 - Compact design so the unit fits on any laboratory table,
 - Very quiet operation,
 - Low unit height for easy insertion and removal of samples,
 - SOFT braking function for slow starting and braking,
 - Optimised for cell separation by gradient centrifugation,
 - "At set RPM" function starts the timer when the selected RPM is reached; to ensure repeatable centrifugation,
 - Push-button lock prevents unintentional parameter adjustment,
 - Stainless steel rotor chamber is easy to clean,
 - Electronic unbalance recognition for maximum safety,
 - Automatic lid lock,
 - Time control: 0.5 - 99 min. and continuous centrifugation,
- Special solutions for centrifuge 5702 R:
 - Temperature range: -9 °C to 40 °C,
 - FastTemp function for rapid cooling,
 - Continuous cooling ensures that the temperature remains constant after operation - samples are cooled all the time,
 - ECO function that switches off the unit after 8 hours of inactivity to reduce energy consumption and extend the life of the compressor,
 - Two program buttons for saving routine programmes,
- Possibility of retrofitting the centrifuge with selected accessories (cups/adapters/samplers) at an additional cost,
- Technical specifications:
 - Max RCF: 3 000 × g,
 - Max speed: 4400 rpm,
 - Max capacity: 4 x 100 ml,
- Included:
 - Pivoting rotor A-4-38B - standard rotor for centrifugation of blood tubes or 15/50 ml conical tubes,
 - Cup for RL120 tubes - 4 pcs,
 - Aerosol - lids for safe centrifugation of hazardous samples - 4 pcs,
 - Adapter 4 x 7 ml - 4 pcs,
 - Adapter 4 x 10 ml - 4 pcs,
 - User manual in Polish (PDF),
 - Power cord,
- Dimensions: 38 x 58 x 26/60 (with open lid) cm,
- Weight: 36 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

