



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **TERMAKS B8054 Laboratory incubator (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Norway,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 330 W,
- TERMAKS B 8054 is a compact incubator that takes up little space with a very high capacity. The interior and shelves are made of stainless steel, The exterior is made of galvanized steel coated with a special gray epoxy paint,
- Specifications:
  - PID temperature control,
  - Max Temp 70°C,
  - Capacity 54 l,
  - Digital display,
  - Audible alarm,
  - Timer,
  - Minute/hour timer: 0-999,
  - Glass door,
  - 4-stage fan speed for fast and better temperature uniformity inside,
  - Door lock,
  - 3 cm access port,
  - Automatic safety thermostat,
- Keypad on control panel in English,
- Includes:
  - 2 x tray,
  - Power cord,
- Dimensions:
  - Overall: 60 x 56 x 57 cm,
  - Compartments: 40 x 33 x 36 cm,
- Weight: 41 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

