



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDIMA S1 Infusion (syringe) pump (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Poland,
- Power supply: 230 V,
- Emergency power supply: through 1 new built-in battery: 9.6 V / 2 Ah,
- Frequency: 50/60 Hz,
- Power: 20 VA,
- Medima S1 is a single-channel syringe infusion pump for precise fluid dosing,
- The pump features:
 - Automatic syringe recognition and grasping,
 - Intuitive operation,
 - Visual and audible alarms with adjustable volume and type of audible signal,
 - Backlit, easy-to-read graphic display,
 - Menu in Polish,
 - Foiled control panel in English,
 - Stand mounting bracket,
 - Wide selection of units,
 - Night mode; reduced display brightness; programmable sound volume, programmable start and end time,
 - Simple connection and disconnection from docking station (without interrupting infusion),
- Infusion types:
 - Continuous,
 - Intermittent,
 - Profile,
- Supported syringes:
 - 5 mL - B. BRAUN OMNIFIX, BD PLASTIPAK, CODAN, ERG Klobuck, KD-JECT III, NIPRO EUROPE, ONCE, TERUMO,
 - 10 mL - B. BRAUN OMNIFIX, BD PLASTIPAK, CODAN, ERG Klobuck, FRES INJECTOMAT, KD-JECT III, MONOJECT, NIPRO EUROPE, ONCE, TERUMO,
 - 20 mL - B. BRAUN OMNIFIX, B. BRAUN PERFUSOR, BD PLASTIPAK, CODAN, ERG KŁOBUCK, KD-JECT III, MONOJECT, NIPRO EUROPE, ONCE, POLFA LUBLIN, TERUMO,
 - 30 mL - B. BRAUN OMNIFIX, BD PLASTIPAK, CODAN, MONOJECT, NIPRO EUROPE, ONCE, TERUMO,
 - 50 mL - B. BRAUN OMNIFIX, B. BRAUN PERFUSOR, BD PERFUSION, BD PLASTIPAK, CODAN, DIPRIVAN, ERG KŁOBUCK, FRES INJECTOMAT, JANPOL, JIANGSU, KD-JECT III, MARGOMED, MONOJECT, NIPRO EUROPE, ONCE, POLFA LUBLIN, SHANCHUAN, TERUMO,
- Specifications:
 - Infusion rate: 0.1 - 2000 mL / h,
 - Infusion volume: 0.1 - 10000 mL,
 - Infusion time: 1 s - 99 h,
 - Bolus dose: 0.1 - 50 mL / h,
 - Bolus infusion: 2 - 2000 mL / h,
 - KVO: 0 - 5 mL / h,
- Dimensions: 36.5 x 23 x 12 cm,
- Weight: 2 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

