



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDIMA DUO HOLDER TP2 Pump docking station (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- DUO HOLDER TP2 is a unique solution for combining single-syringe pumps into two-syringe sets, powered by a single power cord and attached to the stands with a single holder,
- Advantages:
 - Ability to quickly disconnect one of the pumps, e.g., in case it is necessary to use it on another patient,
 - Possibility to create different sets of pumps, including mixed pump sets with peristaltic pumps,
 - Improved infusion safety; track confusion is much less likely than with standard two-syringe pumps; each infusion module is a separate pump, so separate display and keypad,
 - Simple and convenient syringe insertion from the front of the pump,
 - Simple and lightweight design for quick attachment to the tripod,
- The pumps shown in the pictures are not included,
- Overall dimensions: 6.5 x 21.5 x 21.5 cm,
- Weight: 0.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

