



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RODENSTOCK R3 Vision tester (Reconditioned)

- Used,
- Technical condition: good,
- Visual condition: good,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Bulb power: 15 W,
- Optotypes - numbers and pictures,
- Infinitely adjustable head height in the range: 18-30 cm,
- The set includes:
 - Replacement bulb,
 - Test card block,
 - Allen wrench for changing optotype slides,
- Dimensions:
 - Tester: 37 x 18 x 28 cm,
 - Case: 22 x 39 x 34 cm,
- Weight of the set: 7 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украинск

