



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **STRYKER 250 Light source (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- STRYKER 250 is a light source designed to illuminate hard-to-reach areas during procedures, operations or diagnostic procedures,
- Light source: 2 x incandescent bulb 24V / 250 W,
- Features two independent light paths that allow the user to quickly change the bulb during the procedure without interrupting the procedure,
- Lever switch between bulbs for quick change of bulbs,
- Light switch on the front of the device,
- Indicator lights to show the operation of the selected bulb,
- Quiet and efficient cooling system,
- Included:
  - Power cord,
- Dimensions: 41 x 28 x 14 cm,
- Weight: 8 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

