



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **GE 259A KTG unit (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- GE Corometrics 250 series (GE 259A) is a full-featured maternal/fetal monitoring system. It provides comprehensive perinatal monitoring needed during labor. Parameters include the ability to monitor twins, fetal ECG, TOCO/IUP, maternal NIBP, maternal SpO2 and maternal ECG,
- All the best for the baby:
  - Fetal monitoring capabilities include a fetal ECG for displaying and printing continuous, direct measurements of fetal heart rate (FHR) with a curve,
  - Alarms provide both an audible and visual warning when the FHR is outside a user-defined range, in situations where the upper or lower HR threshold is exceeded,
  - An optional and unique feature, Spectra Alerts™, simultaneously analyzes spectral information of FHR and uterine activity (UA) waveforms and alerts physicians to deviations from normal,
- Mother under constant control:
  - Corometrics 250 series provides a choice of integrated maternal pulse oximetry (MSpO2) from companies: GE Ohmeda, Masimo or Nellcor, making it easy for the physician to monitor, document and view the curve,
  - Selectable from a 3-lead maternal ECG signal, with cardiac pulse detection and elimination, it generates the maternal heart rate, displays the QRS curve and a six-second printout segment on demand,
  - Patented DINAMAP SuperSTAT non-invasive blood pressure (NIBP) measurement technology,
- A unique feature introduced by GE, Smart BP delays blood pressure measurement during labor contractions to avoid distracting the patient and obtain more meaningful diagnostic data,
- Nautilus' nine-element waterproof crystalline ultrasound transducers cover a wider focusing region with more uniform coverage at greater depths than conventional seven-crystal transducers. The nine-crystal design allows for minimal transducer repositioning, even for challenging patients of varying body mass and gestational age,
- When monitoring twins, independent ultrasound intensity controls allow the transducer to be applied in such a way that both heart rates can be heard at the same time,
- ChartLight function with new backlit keys with dimmable display,
- Selectable menu in EN, ES, DE, FR, JP, IT, SW, NL,
- Keyboard on the control panel in English,
- Included:
  - User manual in English (PDF),
  - Fetal spiral electrode cable for coronometry (2.4 m),
  - GE Nautilus Toco head,
  - GE Nautilus US head,
  - NIBP cable + cuff,
  - ECG 3-lead cable,
  - ECG bus cable,
  - SPO<sub>2</sub> - brand new silicone sensor,
  - Power cord,
- Dimensions: 45 x 42.5 x 17 cm,
- Weight: 12 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

