



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**TEKNOVA TH-5500 Ultrasound scanner (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in PRC, (2013),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 500 VA,
- TEKNOVA TH-5500 is a high-end fully digital ultrasound imaging system with color Doppler. Equipped with advanced precision DBF (digital beamformer), DRF (dynamic receive focus), RDA (real-time dynamic aperture) and dynamic aperture, the system provides authentic images with rich detail. With a powerful standard PC platform and ergonomic control panel, the system can easily meet the growing demands of clinical diagnostics. With support for a variety of probe types, including linear and convex types, the TH-5500 can be used in many ultrasound diagnostic applications, such as abdominal, gynecology, obstetrics, pediatrics, vascular, cardiology and urology. In addition to excellent image quality, the TH-5500 provides users with various measurement functions. In addition to length, area, perimeter and volume, it also provides special measurements for obstetrics, gynecology, cardiology and urology. The system is equipped with Cine Loop, patient reporting and image management functions. It also supports DICOM 3.0,
- Specifications:
  - 17" high resolution non-interlaced LCD monitor specifically designed for medical imaging,
  - Continuous precision DBF,
  - PC platform,
  - Three active probe ports,
  - Cine Loop min. 512 frames,
  - THI (harmonic tissue imaging),
  - Triplex,
  - SRI (speckle reduction imaging),
  - One-click operating system restore,
  - 320 GB hard drive for permanent storage,
  - DVD-R/W,
  - Two USB ports,
  - DICOM 3.0,
  - Panoramic zoom,
  - Software packages for measurements and calculations,
  - SVVR (high volume video recording for up to 1 hour),
- Imaging:
  - B,
  - B/B,
  - B+M,
  - M,
  - CFM,
  - PW,
  - CDE,
  - CCD,
  - HPRF,
- Selectable menu in the following languages: PL, EN, FR, JP, ES, CN
- Keypad on control panel in English,
- Mobile base on 4 casters with brake,
- Included:
  - Convex C3501 3.5 MHz head,
  - Linear head L7501 7.5 MHz,
  - Micro Convex 3501 3.5 Mhz head,
  - Mitsubishi P93 printer,
  - Power cord,
- Dimensions: 56 x 86 x 140 cm,
- Weight: 72 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN DUTCH  
Мы говорим по русски  
Мы говорим по украински

