



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Toshiba Aplio MX SSA-780A Ultrasound scanner (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1500 VA,
- The Toshiba Aplio MX SSA-780A is a convenient, wide-profile system designed for versatile clinical applications to meet the needs of high-volume medical facilities,
- Specifications:
 - LCD monitor,
 - Screen size: 19",
 - Additional 10.4" touchscreen display,
 - Pull-out alphanumeric keyboard,
 - DVD station,
 - Built-in speakers,
 - 3 active head slots,
 - Transport handle,
 - Maximum imaging depth 28 cm,
 - Frequency range of connectable probes: 1.8 - 18 Mhz,
 - DICOM network,
- Imaging modes: B, M, D, 2D/3D/4D, CW, PW,
- Available features:
 - Pulse Subtraction THI,
 - Differential THI,
 - ApliPure,
 - ApliPure+,
 - QuickScan,
 - MicroPure,
 - Elast-Q,
 - Precision Imaging,
 - Advanced Dynamic Flow (ADF),
- Menu in English,
- Mobile on 4 wheels, all with brake,
- Included:
 - PLT-805AT linear head,
 - Convex head TC375TB,
 - Mitsubishi P95 videoprinter,
 - Power cable,
- Overall dimensions: 85 x 56 x 152 cm,
- Weight: 132 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

