



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **MIDMARK RITTER 230-001 Gynecological chair (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- RITTER 230-001 is a treatment chair designed to make patients feel comfortable and safe,
- Seamless upholstery with medical approval in navy blue, will be replaced with a new upholstery in the color of the customer's choice from our color chart,
- The chair, thanks to its low height and multiple positions, helps to position patients more easily at the touch of a button, reduce the likelihood of fatigue, back and neck strain,
- Infinitely adjustable electrically via foot switch:
  - Footrest angle from 0° to 90°,
  - Seat tilt angle from 0° to 30°,
  - Backrest tilt angle from 0° to 85°,
  - Height of the chair (seat) in the range: 62 - 105 cm,
  - Trendelenburg position: whether to better position the patient, treat shock or fainting, or speed up the onset of anesthesia, the Trendelenburg position provides a full 30° tilt angle,
- Features:
  - Manually adjustable headrest,
  - Extendable, foldable and side-adjustable stirrups,
  - Pull-out treatment bowl,
  - Paper towel storage compartment with width: 45 cm,
- Includes:
  - User manual in English (PDF),
  - Brand new foot switch,
  - Power cord,
- Dimensions:
  - Headrest: 35.5 x 21 x 10 cm,
  - Backrest: 68 x 58 x 10 cm,
  - Seat: 68 x 52 x 10 cm,
  - Footrest: 42 x 37 x 10 cm,
  - Base: 52 x 100 cm,
  - Max. stirrup spacing: 114 cm,
  - Overall length of the chair in reclining position: 185 cm,
  - Overall: 116 x 74 cm x 145 cm,
- Maximum load capacity: 204 kg,
- Weight: 195 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

