



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STIEGELMEYER 6074-03-85-00-0 Bedside cabinet (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Made in Germany,
- Made of particle board finished with veneer,
- Stylish bedside cabinet that can be used in medical centers, hospitals, etc. places,
- Ideal for rooms with limited space,
- Wooden components made of laminated board,
- Double-sided opening drawers and doors that allow the cabinet to be placed on both sides of the bed,
- Easy to clean,
- Rounded edges on all sides of the cabinet, Protects against damage,
- Features:
 - Top with rounded edges with dimensions: 49.5 x 49.5 cm,
 - Drawer equipped with a pull-out separator with dimensions: 40 x 37.5 x 8.5 cm,
 - Open shelf with dimensions: 44 x 40 x 20.5 cm,
 - Lockable space divided into two levels - inside a sliding shelf, allowing the patient to customize,
 - Basket on the door,
- Mobile cabinet, on four plastic wheels,
- Overall dimensions: 50 x 50 x 90 cm,
- Weight: 31 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

