



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE DINAMAP PRO 300 Patient monitor (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Real photos of the product,
- Made in USA,
- Power supply: 230 V and through one new 12V / 2.3 Ah battery,
- Frequency: 50/60 Hz,
- Dinamap PRO 300 - provides non-invasive measurement of blood pressure, heart rate, temperature and oxygen saturation level,
- It has easy-to-use functions and an intuitive menu structure,
- It is extremely versatile and can be used in a wide variety of situations, including acute care facilities such as outpatient clinics, emergency rooms, labor and delivery, and medical/surgical wards,
- The monitor features GE's user interface - making it easy to access patient parameters at the touch of a button,
- Can be used on adult, pediatric and neonatal patients,
- Equipped with:
 - Large easy-to-read liquid crystal display and smaller digital displays,
 - Selectable menu in the following languages: PL, IT, FR, EN, ES, DE, RU, PT, HU, JA, ZH,
 - Built-in printer with paper width: 60 mm,
 - Carrying handle,
 - Trim knob for easy use,
- Product features:
 - DINAMAP NIBP technology - available on GE Solar and Dash monitors,
 - Quality and reliability ix Sigma - built-in quality means accurate readings,
 - Gold Standard Dinamap technology for blood pressure measurement,
 - Flexible mounting options - portable and on wheels,
 - Memory stores up to 100 readings in 24 hours,
 - Protection against water ingress,
 - Intra-arterial reference standard for blood pressure using oscillation technique,
- Specifications:
 - Cuff pressure range:
 - 0 - 290 mmHg (adult),
 - 0 - 140 mmHg (neonate),
 - Heart rate range: 30 - 200 beats per minute,
 - SpO₂ measurement range: 1 - 100%,
- Included:
 - Brand new cuff with air tubing,
 - Brand new SpO₂ clip sensor from BIOGENESIS,
 - User manual in English (PDF),
 - Power cord,
- Dimensions:
 - Overall: 49 x 49 x 122 cm,
 - Monitor: 18 x 26 x 27 cm,
- Weight of the set: 14.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

