



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SIEMENS CLINITEK STATUS+ Urine analyzer (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Realistic pictures of the product,
- Made in USA,
- Power supply: 9 V or battery operated (6 x AA/1.5V),
- Frequency: 50/60 Hz,
- Power: 7.2 VA,
- Clinitek Status+ is an analyzer whose most important feature is a system equipped with Auto-Checks*. The unit automatically checks each test strip for moisture exposure, typical interferences found in samples, and recognizes the Siemens strip type. All this together improves the quality of the clinical information obtained,
- Features:
 - Reporting results in clinical, conventional or SI units,
 - Marking pathological results,
 - Detection of fresh and leached blood cells,
 - Ability to determine urine color and clarity,
 - Free display of results stored in memory,
 - Comfortable touch screen,
 - Selectable menu in the following languages: EN, DE, FR, ES, IT, SV, ZH, JA,
 - Choice of two test modes: rapid test and extended test,
 - Ability to enter patient information, identification, operator and password to access the instrument settings,
 - Ability to transfer data to a computer (RS232 socket),
 - Test panel expanded to include rapid immunoassay,
 - Has a built-in thermal printer,
- Supported strips:
 - Hema-Combistix,
 - Uristix,
 - Uristix 4,
 - N-Neostix,
 - Uro-Uristix,
 - Combistix,
 - Multistix 10 SG,
 - Multistix 8 SG,
 - Multistix 7,
 - Multistix 5,
 - N-multistix SG L,
 - Uro-Hema-Combistix SG-L,
 - Uro-Hema-Combistix,
 - Labstix SG,
 - Multistix SG,
 - N-Multistix SG,
 - Multistix GP,
- Specifications:
 - Performance: approximately 60 studies per hour,
 - Memory: 700 test results including patient data,
 - Automatic calibration, performed automatically,
- Included:
 - User manual in Polish and English (PDF),
 - Power cord with power adapter,
- Dimensions: 16.5 x 27 x 16 cm,
- Weight of the set: 2 kg,
- Cat 01,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN DUTCH
Мы говорим по русски
Мы говорим по украински

