



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE LOGIQ P5 Ultrasound scanner (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Real photos of the product,
- Made in USA, (2011),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 750 VA,
- GE Logiq P5 is a powerful ultrasound machine with shared services. It supports a variety of applications, including OB/GYN and cardiology, and includes many features not usually found in a mid-range system. It can perform 4D obstetric scans (optional), support surgical procedures and perform a variety of cardiac scans. Premium features of the GE Logiq P5 include SRI (speckle reduction imaging), CE (coded excitation), ASO (automatic spectrum optimization), Crossbeam technology and more. The GE Logiq P5 is designed with a simple, ergonomic layout that enables fast workflow and efficiency. It features a 15-inch TFT monitor with a fully movable arm for better viewing angles for sonographers and patients,
- Specifications:
 - 15-inch LCD monitor on an articulating arm,
 - Full-size keyboard with user-defined keys,
 - Automatic Optimization (AO),
 - Speckle Reduction Imaging (SRI),
 - 3D/4D imaging (appropriate head required),
 - TrueScan architecture,
 - True Access raw data processing,
 - SmartScan,
 - ComfortScan,
 - Integrated hard drive and DVD-R/W,
 - Built-in memory for peripheral devices,
- Imaging modes:
 - B-Mode,
 - Harmonic Imaging,
 - M-Mode,
 - Color Flow Mode (CFM),
 - Power Doppler Imaging (PDI) ,
 - PW Doppler,
- Applications:
 - Obstetrician/Gynecologist,
 - Vascular,
 - Abdominal,
 - CEC,
 - Urology,
 - Cardiology,
 - Pediatrics,
 - Neonatal,
 - Transcranial,
- Selectable menus in the following languages: UK, DE, IT, FR, ES, RU, DK, NL, CN, JP,
- Keypad on control panel in language: English,
- Mobile on 4 wheels with brake,
- Included:
 - User manual in English (PDF),
 - Convex head 4C,
 - Linear head 12L,
 - SONY UP-D897 Videoprinter,
 - Power cord,
- Overall dimensions: 52 x 62 x 136 cm,
- Weight: 80 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN DUTCH
Мы говорим по русски
Мы говорим по украински

