



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG Type VACUKLAV 31 Autoclave (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2500 W,
- MELAG Vacuklav 31 is an efficient steam autoclave offering 4 work programs and 2 test programs,
- Product features:
 - Class B,
 - Automatic,
 - Pre-vacuum and fractional vacuum produced by vacuum pump,
 - Possibility of additional vacuum drying,
 - Preheat option helps maintain temperature between cycles,
 - Equipped with indicator lights to inform you of unit errors,
 - 2.5" LCD display,
 - Foiled control panel,
 - Keypad on the control panel in German,
 - Place for a tray on the lid of the autoclave,
 - Water drain connector on the front of the unit,
 - Sockets:
 - For Melaprint printer (no printer included),
 - For Meladem ion exchange system (system not included),
 - 230 V outlet,
 - Made of stainless materials,
 - Stable on legs, including front legs with height adjustment,
- Features 4 work programs and 2 test programs:
 - Unpacked tools: pressure 2 bar - time 4 min - temperature 134 °C,
 - Packaged tools: pressure 2 bar - time 5 min - temperature 134 °C,
 - Textiles: pressure 1 bar - time 17 min - temperature 121 °C,
 - Disinfection: pressure 0.3 bar - time 8 min - temperature 105 °C,
 - Bowie & Dick test: vapor permeation test to check the efficiency of the autoclave by the user,
 - Leakage test: vacuum test to check the leakage of the autoclave,
- Specifications:
 - Sterilization chamber capacity: 14,5 L,
 - Maximum load capacity: 4 kg of instruments including trays or 1.5 kg of textiles,
 - Water consumption during sterilization: 300 - 400 mL, depending on the program and loading level,
- Included:
 - Sieve trays with dimensions: 18 x 28 x 2 cm - 2 pieces,
 - Tray holder,
 - Stand with holders for 3 trays,
 - Brand new condensate tank,
 - User manual in English (PDF),
- Dimensions:
 - Autoclave: 63 x 45 x 45 cm,
 - Chamber: Ø 23 cm x 34 cm,
- Weight: 44 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN DUTCH
Мы говорим по русски
Мы говорим по украински

