



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Bedside cabinet (Reconditioned)

- Used,
- Technical condition: good,
- Visual condition: good,
- Bedside cabinet is an ideal solution for rooms with limited space,
- Product features:
 - The painted cabinet body guarantees high resistance to wear, corrosion and peeling,
 - The work surface and eating boards are made of high-pressure laminate with high mechanical, chemical and thermal resistance,
 - The eating board is equipped with raised edges to prevent the object from falling,
 - Smooth surfaces for easy cleaning,
- Equipped with:
 - Extendable and tiltable side table top with adjustable height from 83 - 119 cm, which allows the patient to position it at different angles for activities such as eating, reading,
 - A small top drawer with pull-out plastic waterproof organizers,
 - An alcove for storage,
 - A lower cabinet that opens from one side,
 - Shoe shelf under the cabinet,
 - Four double casters (including two with brakes),
- Dimensions:
 - Working surface of the upper top: 45 x 46 cm,
 - Eating board: 42 x 52 cm,
 - Top drawer: 29 x 32.5 x 6.5 cm,
 - Recess: 36 x 39.5 x 8 cm,
 - Bottom drawer: 34 x 39 x 28 cm,
 - Overall: 45 x 72 x 80 cm,
- Weight: 27.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

