



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CHIRANA MEDICAL SK1 Dental chair with lamp (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Slovakia (2007),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 300 VA,
- **Chirana SK1** is a compact, office space-saving dental chair with a lamp,
- **New medically certified upholstery in gray (upholstery can be replaced, from our color pattern book for additional fee),**
- **Smooth, electric regulation with the use of a footswitch:**
 - Seat height in the range of 45 - 92 cm,
 - Backrest angle in the range of 0° / 70°,
- **Smooth, manual regulation:**
 - Headrest angle,
 - Headrest height in the range of 8 cm,
- **The chair is equipped with:**
 - **CHIRANA dental lamp:**
 - Mounted on a column with a pantograph arm,
 - Light source: halogen bulb 17 V / 95 W,
 - Two-stage regulation of light intensity,
 - Manual multi-level regulation of the lamp setting, including the height of the position of the canopy within 40 cm,
 - Anatomical headrest for maximum patient comfort,
 - Left armrest,
- **In the set:**
 - Two-functional footswitch in the form of a joystick,
- Dimensions:
 - Of the chair: 150 x 48 x 121 cm,
 - Overall: 150 x 65 x 196/241 cm,
- Maximum load: 200 kg,
- Weight: 134 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

