



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



PENLON PRIMA SP2 Anesthetic machine (Respirator/Ventilator) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in UK,
- Power supply: 230 V,
- Emergency power supply: via 1 new built-in 12 V / 1.3 Ah battery,
- Frequency: 50/60 Hz,
- Power: 1000 - 2500 VA,
- The Penlon Prima SP2 is a machine for anaesthetising the patient during surgery,
- Features two ports for vaporisers with common Dräger connector systems (no vaporisers included - can be retrofitted at additional cost),
- Supplied with gases from a central system: O₂, N₂O, AIR,
- Can be supplied with gases from additional O₂ or N₂O cylinders,
- Respirator driven by an electric motor (electronically controlled),
- Compact patient breathing system with oxygen sensor located on a mobile arm,
- Has a monitor on the mobile arm,
- Selectable menu in the following languages: PL, EN, CE, ES, FR, IT, PT,
- Apparatus equipped with:
 - Work surface with dimensions: 62.5 x 27 cm,
 - Three drawers with dimensions: 52 x 35.5 x 15 cm,
 - A rail for attaching additional accessories,
 - Light for illuminating the worktop,
- For an additional charge, the unit can be retrofitted with:
 - Vital function monitor / cardiomonitor,
 - Injector suction unit,
- Mobile apparatus, on 4 double castors, 2 of them with brake,
- Included:
 - Gauze supply hoses with tips of the customer's choice (AGA or DIN),
 - Patient breathing system,
 - User manual in English (PDF),
- Overall dimensions: 70 x 110 x 136 cm,
- Maximum load capacity of drawers: 10 kg,
- Weight: 133 kg,
- In the domestic market (Poland), in the event of installation of the equipment within 30 days of delivery, installation and training are free of charge, after this deadline these services will be charged,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

