



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS OEV191 LCD medical monitor (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 125 VA,
- The Olympus OEV191 is a large medical monitor delivering high quality images,
- Progressive scanning reduces image distortion and flicker,
- Anti-glare coating,
- Image enhancement function,
- 19" TFT LCD matrix,
- Horizontal and vertical viewing angles: 85°,
- Resolution: 1280 x 1024 px,
- Menu language: English,
- Keypad on control panel in English,
- The monitor has slight signs of use,
- Features following inputs and outputs:
 - VIDEO IN-OUT, Y/C IN-OUT,
 - RGB/COMPONENT: A (G/Y IN-OUT, B/Pb IN-OUT, R/Pr IN-OUT, EXT SYNC IN-OUT, B (G/Y IN, B/Pb IN, R/Pr IN, EXT SYNC IN)
 - REMOTE: RS - 232, REMOTE 1,
- Dimensions: 43 x 39 x 10 cm,
- Weight: 8 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 1 month for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински



