



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ALOKA PROSOUND SSD-3500SX Ultrasound scanner (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Japan,
- Power supply: 320 V,
- Frequency: 50/60 Hz,
- Power: 700 VA,
- The Aloka Prosound SSD-3500SX offers clear imaging, high versatility, user-friendly design and comprehensive data management,
- Specifications:
 - 15" LCD screen,
 - 3 active head ports,
 - Transport handles,
 - Alphanumeric keyboard built into the control panel,
 - Built-in speakers,
 - USB 2.0 port,
 - Advanced measurement/reporting features,
- Imaging modes:
 - M,
 - B,
 - D,
 - B/B,
 - B/M,
 - B/D,
 - Harmonic,
 - Flow,
 - Power Flow,
 - 2D,
 - Color Doppler,
 - Power Doppler,
- Menu in English,
- Keypad on control panel in English,
- Mobile on 4 wheels, all with brake,
- Included:
 - Aloka convex head UST-9123,
 - Aloka endovaginal head UST-9124,
 - Printer Mitsubishi P93,
 - Power cord,
- Overall dimensions: 49 x 79 x 149 cm,
- Weight: 128 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

