



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **FAMED ŁÓDŹ LU301E Negatoscope (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- Real photographs of the product,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 120 VA,
- Famed LU301E is a device for viewing X-rays, slides, drawings on film, tracing paper, etc.,
- Durable metal casing with easy disinfection of the equipment, has visible signs of use,
- Easy servicing (replacement of light sources),
- Even illumination of the screen,
- Advantages:
  - Bright light close to daylight, for viewing X-ray and CT images,
  - Infinitely adjustable screen illumination intensity by means of a knob,
  - Very easy replacement of worn out fluorescent tubes and ignition system components,
  - Convenient and reliable fixing of film of different thicknesses and dimensions,
- Specification:
  - Light source: 3 x 36 W/ 54 (120 cm),
  - 1-section, without individual frame illumination,
  - Negatoscope can be hung on the wall,
- Dimensions:
  - Overall: 13 x 47 x 134.2 cm,
  - Screen: 40 x 118 cm,
- Weight: 14.5 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

