



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **OLYMPUS EVIS EXERA III CF-H185L Video colonoscope (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- The CF-H185L with EVIS EXERA III provides HDTV-quality images for more detailed observation of mucous membranes and capillaries,
- Advantages:
  - NBI technology allows observation from twice the distance of the EVIS EXERA II 180 series endoscopes, while providing higher contrast between blood vessels and mucosa,
  - Adjustable stiffness avoids retraction of the endoscope, depending on the specific case, to accommodate the patient's unique anatomy and the physician's preference for manipulating the instrument,
  - The zoom function provides a magnified image of the selected point by positioning the endoscope tip as close as 2 mm from the mucous membranes,
- Specifications:
  - Field of view: 140°,
  - Viewing direction: forward view,
  - Depth of field: 2 - 100 mm,
  - Outer diameter of distal tip: 12.8 mm,
  - Outer diameter of insertion tube: 12.8 mm,
  - Working length: 1680 mm,
  - Internal canal diameter: 3.7 mm,
  - Minimum visible distance: 3.0 mm from the distal end,
  - Tip bending angles:
    - Upwards: 180°,
    - Downwards: 180°,
    - Right: 160°,
    - Left: 160°,
  - Overall length: 2005 mm,
  - Compatible with EVIS EXERA III (190) systems,
- Included:
  - Valves,
  - Adapter for suction and washing of canals (MH-856),
  - Adapter for manual canal decontamination (MH-946),
  - Biopsy channel inlet and valve seats plug (MH-944),
  - User manual in English (PDF),
- Case dimensions (carton): 60.5 x 42 x 13 cm,
- Weight of set: 4 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински



PRAISTON

