



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **OLYMPUS GIF-H185 Video gastroscope (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan (2015),
- OLYMPUS GIF-H185 is a video gastroscope that provides the best possible HDTV image quality,
- Product features:
  - The EVIS EXERA III system guarantees HDTV-quality images, i.e. high-resolution viewing capability, even with smaller diameter endoscopes,
  - NBI technology (narrow band imaging) - enables observation from twice the distance of the EVIS EXERA II 180 series endoscopes, while providing higher contrast between blood vessels and mucosa,
  - Close-ups - The close-up function allows a magnified image of a selected point by positioning the tip of the endoscope as close as 2 mm from the mucous membranes,
  - Slim design - This endoscope stands out for its excellent balance of size and performance. HDTV image quality in a slim size with a diameter of 9.2 mm,
  - Waterproof single-stage connector - Reduces the list of set-up tasks performed before and between examinations of successive cases. In addition, it can be completely submerged in water, eliminating the need for a waterproof cap and the risk of costly repairs due to accidental immersion,
- Technical specifications:
  - Optical system:
    - Field of view: 140°,
    - Viewing direction: forward view,
    - Depth of field: 2 - 100 mm,
  - Endoscopic probe:
    - Outer diameter of endoscopic probe tip: 9.2 mm,
    - Outer diameter of endoscopic probe tube: 9.2 mm,
    - Working length: 1030 mm,
  - Biopsy channel:
    - Inside diameter of channel: 2.8 mm,
    - Minimum viewing distance: 3.0 mm from the tip of the endoscopic probe,
  - Tip bendable in the range:
    - 210° upwards,
    - 90° downwards,
    - 100° clockwise,
    - 100° left,
- Compatible with EVIS EXERA system:
  - OLYMPUS CV - 190 video processor,
  - OLYMPUS CLV - 190 xenon light source,
- Included:
  - Valves,
  - User manual in English (PDF),
- Overall length: 135 cm,
- Case dimensions (cardboard box): 60 x 42 x 13 cm,
- Weight of set: 4 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

