



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERON Medical couch (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- Real photographs of the product,
- Made in Germany,
- ERON is a two-segment couch designed for treatment rooms, doctors' surgeries, specialist and therapy rooms for medical examinations and minor medical procedures, massage and rehabilitation,
- Upholstery with medical certificate in grey colour - has traces of use visible in the last photos (possibility of replacing the upholstery with a new one from our colour chart for an additional charge),
- Sturdy construction, powder coated in white,
- Specifications:
 - Wide two-section recliner,
 - Construction made of powder-coated steel profiles,
 - Manual stepped headrest adjustment with a knob in the range: 0 - 56°,
 - Has a paper hanger with a width of: 70 cm,
- Overall dimensions: 197 x 82.5 x 66 cm,
- Weight: 42.5 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

