



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LIFEPAK CR2 Training defibrillator (Ex-Display)**

- New (ex-display),
- Technical condition: perfect,
- Visual condition: very good,
- Real photos of the product,
- Made in USA, (2018),
- Power supply: via four 1.5 V alkaline D (R20) batteries,
- The LIFEPAK CR2 training defibrillator enables training under realistic conditions without the need to charge and discharge the device with electricity. The training version looks exactly like the clinical LIFEPAK CR2, is controlled by the same simple user interface, features all the functions of a clinical device available, so it perfectly simulates its actual use. Perfect for training with any manikin or phantom model without the need for additional accessories,
- Simple, intuitive to use and affordable, the training equipment fully replicates realistic conditions, guiding users through simulated analysis, defibrillation performance and correct CPR cycles with support messages,
- Like the clinical CR2, the training version uses ClearVoice technology to detect noise and adjust the volume level for intelligible instructions,
- The training defibrillator can be programmed to simulate a variety of scenarios that users may encounter when using the LIFEPAK CR2 defibrillator in real-world conditions. The training device has two completely customisable scenarios pre-programmed, which can be fully adapted to the situation on the settings screen and modified for maximum flexibility,
- With the CR2 Trainer Demo app, the function can be set up to automatically receive email notifications or send manually prepared email messages via Bluetooth from any Windows, Android or IOS mobile device. The app provides a simulation of automatic email notifications sent from CR2 devices that are connected to the LIFELINKcentral AED programme manager,
- It detects when pressure is applied to the QUIK-STEP™ electrodes and automatically cycles through the steps of the scenario, eliminating the need for remote control. A pause button on the device, if required, allows the scenario in progress to be stopped,
- The CR2 AED training system is cost-effective. It is powered by four standard D-cell batteries and the electrodes are affordable and reusable and are an exact replica of the QUIK-STEP defibrillation ECG electrodes available quickly and used exclusively on the LIFEPAK CR2,
- To enable a quick demonstration of the basic functions of the LIFEPAK CR2 Defibrillator and the benefits of cprINSIGHT analysis technology, scenarios can be programmed for automatic or semi-automatic mode,
- Selectable menu in the following languages: PL, DE, EN, ES, SE, CN, JP, DK, NL, FR, IT,
- Included:
  - Electrode set,
  - 4 alkaline D (R20) 1.5 V batteries,
  - Transport bag for defibrillator with strap,
  - User manual in Polish (PDF),
- Dimensions: 23 x 27.5 x 10 cm,
- Weight: 3 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

