



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HILL-ROM AVANTGUARD 1200 Hospital bed (ICU) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in France,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V / 1.2 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 360 VA,
- The AvantGuard 1200 is a multifunctional hospital bed designed to meet the needs of the entire medical team,
- Due to its functions, it is designed for use not only in general care, but also in post-operative care and in the ICU,
- It facilitates the use of vital signs monitoring devices and patient transport,
- Has buttons for selective function locking located on both sides of the bed,
- Comes with a new mattress shown in the last photo,
- Electrical infinite adjustment via the control panel:
 - Bed height in the range: 55 - 95 cm,
 - Back section angle in the range: 0 - 65°,
 - Trendelenburg and anti-Trendelenburg positions: + 16° / - 16°,
 - Tilt angle of thigh segment: 36°,
 - Foot segment tilt angle: - 23°,
- Equipped with:
 - Safety folding handrails with holes,
 - Central brake,
 - 4 drip spaces in each corner of the bed,
 - 4 corner bumpers,
 - Double-sided levers for CPR,
 - Hand controls for the carer,
 - Removable head and leg tops,
 - 8 attachments for patient straps,
 - Integrated bed extension,
 - Bed linen shelf,
 - Brake lock alarm,
 - Pedals for height adjustment on both sides,
- X-RAY translucent,
- Mobile, on 4 double wheels with Ø 15 cm,
- Has normal signs of use,
- Included:
 - Mattress with dimensions: 195 x 85 x 13 cm,
 - User manual in English (PDF),
- Overall dimensions: 212 x 99 x 96 cm,
- Load capacity:
 - Safe: 185 kg,
 - Maximum: 250 kg,
- Weight: 185 kg,
- Cat 03,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

