



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARL ZEISS OPMI-9 surgical microscope (ENT) (wall-mounted) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 70 VA,
The OpMi-9 is a compact, easy-to-use microscope,
Illumination - bulb: 6 V / 30 W,
Optical system:
Eyepiece with magnification: 12.5x with adjustable focus,
Objective with focal length: $f = 250$ and $f = 500$,
Manually adjustable interpupillary distance,
Key design features:
Device suitable for wall mounting,
Three-jointed arm with height adjustment in the range of: 35 cm,
Rotating arm with microscope body,
Adjustable in many planes,
Maximum span of the set: 90 cm,
Dimensions Overall dimensions: 30 x 100 x 90 cm,
Set weight: 20 kg,
Valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
3 months for the international market,

