



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **CARL ZEISS OPMI-9 surgical microscope (ENT) (wall-mounted) (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Germany,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 70 VA,  
The OpMi-9 is a compact, easy-to-use microscope,  
Illumination - bulb: 6 V / 30 W,  
Optical system:  
Eyepiece with magnification: 12.5x with adjustable focus,  
Objective with focal length:  $f = 250$  and  $f = 500$ ,  
Manually adjustable interpupillary distance,  
Key design features:  
Device suitable for wall mounting,  
Three-jointed arm with height adjustment in the range of: 35 cm,  
Rotating arm with microscope body,  
Adjustable in many planes,  
Maximum span of the set: 90 cm,  
Dimensions Overall dimensions: 30 x 100 x 90 cm,  
Set weight: 20 kg,  
Valid Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
3 months for the international market,

