



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FOTONA STARWALKER M031-3A/2 MAQX Laser (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Slovenia, (2018),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- The FOTONA StarWalker laser is a modern third-generation laser, a combination of up to 4 laser wavelengths and 14 modes of operation on different pulse types for maximum effectiveness. FOTONA StarWalker will effectively deal with even the most difficult tattoo pigmentation, hyperpigmentation and vascular lesions without the risk of damaging the skin surrounding the lesion,
- Possible applications:
 - Removal of scars and acne scars,
 - Elimination of stretch marks,
 - Removal of hyperpigmentation,
 - Treatment of acne and rosacea,
 - Removal of tattoos,
 - Anti-aging treatments,
 - Laser treatment of nail fungus,
 - Vascular removal,
 - Removal of erythema,
- Included:
 - R28 d head (with a round spot size of 2 - 8 mm),
 - Head R58 d (with round spot size 1.5 - 6 mm),
 - FS20 d head (1064 nm 9mm 5x5),
 - FS50B head d (532 nm 5mm 5x5),
 - Foot switch,
 - Power cable,
- Weight: 68 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

