



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ophthalmic table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- Real photos of the product,
- Ophthalmic table with manual height adjustment,
- Specifications:
 - Table top 53 x 34 cm,
 - Height adjustment in the range: 70-90 cm,
 - Base diameter: 45 cm,
 - Top made of laminated board - has normal signs of use including a few small holes from fixing screws,
 - Powder-coated base made entirely of metal,
 - Shelf under the top with internal dimensions 30.5 x 7.5 x 11 cm,
- Overall dimensions: 50 x 53 x 70/90 cm,
- Weight: 32.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

