



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **B.BRAUN PERFUSOR SPACE Infusion (syringe) pump (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V and through a new built-in 1.8 Ah NiMH battery,
- Frequency: 50/60 Hz,
- Power consumption: 8 VA,
- Perfusor Space is a syringe infusion pump for intravenous, intra-arterial infusion and enteral feeding and transfusion,
- Specifications:
  - Operating time when battery powered: about 8 hours of infusion at 25 ml/h rate,
  - Flow rate: 0.01 - 200 ml/h,
  - Bolus: 1800 ml/h,
  - Time preselection: 00:01 - 99:59 hrs,
  - Standby mode: 1 minute ... 24 hours, drug name also visible in standby mode,
  - Drug content - amount of active agent in: µg, mg, g\*, mEq, IU, mmol, kcal, kJ,
  - Visual and acoustic alarms,
  - Barcode assignment function,
- Foiled and backlit keyboard in English,
- Backlit graphic display, visibility angle ~40° on all sides,
- Selectable menu in English and German,
- At additional cost:
  - Possibility to change language to: Hungarian, Swedish, Portuguese, Norwegian, Italian, Dutch, French, Finnish, Spanish, Danish, Czech, Polish, Russian, English (US), Turkish, Bulgarian, Greek, Indonesian, Romanian,
  - Possibility to change the software to support syringes: Monoject 50ml EU, Monoject 60cc US, Codan 50/60ml, Polfa 50ml, Margomed 50ml, Injectomat 50ml,
- The pump supports syringes:
  - BBraun OPS 20,
  - BBraun OPS 50 ml,
  - BBraun OMNIFIX 10 ml,
- Included:
  - Brand new power supply,
  - User manual in English (PDF),
- Dimensions: 26.5 x 20.5 x 10 cm,
- Weight:
  - Pump: 1.5 kg,
  - Set: 1.8 kg,
- Cat 03,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

