



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DATEX-OHMEDA ISOTEC 5 Vaporiser for anaesthetic agents (Isoflurane) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in USA,
- GE Datex-Ohmeda TEC 5 Isoflurane is designed for off-circuit use in continuous flow techniques for inhalational anaesthesia. TEC 5 is temperature, flow and pressure compensated so that its performance remains relatively constant despite cooling due to evaporation,
- The TEC 5 evaporator can be installed in Datex-Ohmeda anaesthesia machines equipped with the Selectatec attachment manifold. No tools are required for this operation and the anaesthesia machine does not have to be taken out of service,
- The liquid anaesthetic level can be read thanks to an indicator on the front of the evaporator,
- Liquid capacity (dry wick): 300 ml,
- Liquid retained by the wick: 75 ml,
- Temperature operating range: 18°C to 35°C,
- Has normal signs of use including a small dent on the control knob (photo) not affecting the vaporiser's work in any way,
- Included:
 - User manual in English (PDF),
- Dimensions: 25 x 11 x 6.5 cm,
- Weight: 7 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

