



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **REHA FUND BIOFLOTE 2000 Anti-decubitus mattress with pump (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in PRC, (2019),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Pneumatic anti-decubitus mattress with bubble design for short-term home care. Ideal for people with risk of pressure sores or pressure sores up to grade I. It is dedicated to short-term home care during post-operative periods, during rehabilitation or for extended periods of lying down for most of the day. The mattress is constructed with 130 mini-chambers and two main sections, which alternately inflate or deflate, reducing the pressure of the base on the entire body surface of the person lying down in a 12-minute cycle. This arrangement allows the blood vessels at the pressure points on the body to be unclogged, thus supplying oxygen and nutrients to the soft tissues with blood. The quiet operation of the pump, the ease of adjusting the pressure in the mattress in relation to the patient's weight body and the sturdy construction are just some of the advantages of this mattress, which is one of the most popular models on the market,
- Specifications:
  - Quiet pump (<31dB),
  - 2 sections alternately inflate or deflate in a 12 min cycle perfectly redistributing pressure,
  - Padding at the ends prevents the mattress from slipping,
  - A convenient gauge on the pump with a range of patient weights, making mattress pressure selection completely automatic,
  - Non-slip hooks on the pump,
  - Hooks for hanging the pump included,
  - Hose made of durable and soft medical-grade PVC (polyvinyl chloride),
  - Easy to clean,
  - Easy to use,
  - Easy to assemble,
  - Cradle made of highly flexible, medical-grade PVC,
- Dimensions:
  - Mattress when inflated: 200 x 90 x 6.5 cm,
  - Pump: 27 x 11 x 9.5 cm,
- Maximum patient weight: 130 kg,
- Weight: 3 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

