



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Certified CLARO MEDTECH Sars-CoV-2 & Influenza A/B Antigen combo test/tests - 19 (1 pc/package) (New)**

- New,
- Real photos of the product,
- Made in PRC,
- Comes with an approved CE mark,
- Sars-CoV-2 & Influenza A/B Antigen is an immunochromatography-based one-step in vitro test,
- It uses a dual antibody sandwich method for the detection of SARS-CoV-2 and influenza A/B by colloidal gold immunochromatography,
- It is designed for the rapid qualitative determination of SARS-CoV-2 antigens, influenza A and B viruses in anterior nasal swabs taken from individuals with suspected COVID-19, influenza A and B within the first seven days of symptom onset,
- The test is intended for persons aged 18 years or older with clinical signs of SARS-CoV-2 and influenza A and B infection or suspected SARS-CoV-2 influenza A and influenza B infection,
- The test kit is for self-testing (layman's test),
- Specifications:
  - Response time: 15 - 20 minutes,
  - Sample: nasal swab,
  - Storage : 2 - 30° C,
- Included:
  - Test cassette in sealed plastic bag (with absorbent) - 1 pc,
  - Lysis buffer and dropper (350 µl / tube) - 1 pc,
  - Disposable sterile packaged anterior nasal swab stick - 1 pc,
  - User manual in English,
  - Biological bag - 1 pc,
  - Certificate,
- Packaging dimensions: 2 x 15 x 5.5 cm,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

