



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CUDA ELITE XLS-300 Light source + CUDA DRH-110 Headlamp (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 450 W,
- The CUDA Elite XLS-300 light source provides the brightest and whitest illumination required for surgical procedures,
- Light source: Xenon 300W SYC0052,
- Colour temperature 5600K,
- Infinitely variable control of light intensity: 0 - 100 %,
- Available ports: ACMI, Wolf, Storz, Olympus,
- Rear panel timer,
- Included:
 - User manual in English (PDF),
 - CUDA DRH-110 headlamp,
 - Power cable,
- Dimensions: 13 x 25 x 36 cm,
- Weight: 7 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

