



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **KAMPLEX AD12 Diagnostic audiometer (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- Real photographs of the product,
- Made in England,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- The AD12 is a portable audiometer that is ideal for schools, industry and primary care practices or wherever rapid hearing screening is carried out,
- Easy to use, requires minimal training prior to use and offers calibrated pure tone and tone stimuli over a wide frequency range,
- Specifications:
  - Frequency range 250 - 8000 Hz,
  - Volume range from -10 to 100 dB,
  - Control panel keypad in English,
- Included:
  - Re-loop RH-2350 PRO patient headphones,
  - Original carrying case,
- Dimensions:
  - Audiometer: 24 x 36 x 12.5 cm,
  - Case: 16 x 39 x 41 cm,
- Weight of set: 4.5 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

