



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MALLINCKRODT PURITAN-BENNETT 840 Ventilator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in USA,
- Power supply: 230 V and rechargeable battery,
- Frequency: 50/60 Hz,
- The Puritan-Bennet 840 is a ventilator for long-term replacement ventilation for intensive care purposes,
- It can meet the specific needs of virtually any type of patient,
- It automatically compensates for the flow-dependent drop in endotracheal tube or tracheal pressure and controls the patient's cerebral pressure to a constant PEEP value,
- When the circuit is disconnected, the gas flow stops and the clinicians are alerted,
- Ongoing background checks assess the ventilator's electronics and ventilation to ensure patient safety,
- Innovative design makes setup and navigation easy,
- The ventilator supports hospital information systems and provides seamless data transfer,
- DualView LCD touchscreen,
- Menu and control panel in English,
- The ventilator is equipped with the following software:
 - BiLevel TM - allows the patient to breathe spontaneously at two levels of PEEP,
 - TC (Tube Compensation) - Compensating for the resistance in the artificial airway, this type of spontaneous breathing exactly overcomes the imposed work of breathing (WOB) from the endotracheal and tracheotomy tube,
- Ventilator parameters:
 - Pressure support (PSUPP): 0 to 70 cm H₂O,
 - Expiratory sensitivity range (ESENS): 1% to 80%;
 - Apnea tidal volume velocity (VT): 25 to 2500 ml,
 - Apnoea respiratory rate (f): 2.0 to 40 / min,
 - Peak inspiratory flow (V.MAX):
 - 1.0 to 30 l / min, for neonates,
 - 3.0 to 60 l / min, for pediatric patient,
 - 3.0 to 150 l / min for an adult patient,
 - Inspiratory pressure (PI): 5 to 90 cm H₂O,
 - Apnoeic inspiratory time (TI): 0.2 to 8.0 seconds,
 - Apnoea interval (TA): 10 to 60 seconds,
 - Apnoeic expiratory time (TE): TE ≥ 0.2 seconds,
 - End-expiratory pressure: PEEP: - 20 to 130 cm H₂O,
- Included:
 - KORGIEL Model W802 Oxygen reducer,
 - Oxygen and air supply hoses with customer's tips choice (AGA or DIN),
 - Adult breathing circuit with 2.0 L breathing bag,
 - Reusable filter with autoclavable condensate container,
 - User manual in Polish and English (PDF),
- Mobile on four castors, two with brakes,
- Overall dimensions: 58.5 x 60.5 x 154 cm,
- Weight: 47.5 kg,
- Cat 08,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

