



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BERCHTOLD CHROMOPHARE D200 Diagnostic and treatment light (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 45 VA,
- BERCHTOLD D200 is a wall-mounted diagnostic and treatment lamp,
- Specifications:
 - Built-in transformer,
 - Sterilizable handle,
 - Light source: 12 V 50W incandescent bulb,
 - Infinitely adjustable light intensity,
 - Max light intensity 25,000 Lux at a distance of 0.8 m,
 - Color temperature: 3600 K,
 - Lamp range: 60 -130 cm,
- Dimensions:
 - Total length of the lamp about 150 cm,
 - Diameter of the canopy: 21 cm,
 - Reflector diameter: 16.5 cm,
- Weight: 4 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

