



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MELAG 75 Dry heat sterilizer (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 690 W,
- MELAG 75 is a dry heat sterilizer, a fully safe and durable device. It meets all the requirements for compact devices,
- Safety is guaranteed by the function of automatic shutdown of the sterilization cycle including power off when the cycle is completed,
- The device can be used in veterinary and medical fields,
- Dry sterilization,
- Forced air circulation,
- Capacity: 4.3 L,
- Infinitely adjustable temperature from 60° C to max. 180° C,
- Has a timer switch adjustable from: 0 - 120 min,
- Stainless steel housing, strong construction and long device life,
- Includes:
  - Tray,
  - Tray holder,
  - User manual in German (PDF),
  - Power cord,
- Dimensions:
  - Tray: 18 x 28 x 2 cm,
  - Chamber: 19.5 x 29 x 7.3 cm,
  - Overall: 31 x 41.5 x 26 cm,
- Maximum load of tray: 2 kg,
- Weight: 12 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

